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BIBLIOGRAFIA

Gorge P. Schmid. Current Opinion in Infectious Diseases 1994;7 34-40.

Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.

Sandra Larsen et al. A manual of test for Syphilis American Public Health Association 190: 1-192.

Joseph Earle Moore et al. Gastrointestinal Haemorrhage 1952; 150 (5): 467-473.


Young DS. Effects of drugs on clinical laboratory test, 4th ed AACC Press. 1995.

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Qualitative determination of plasma reagins

PRINCIPLE

The RPR carbon is a non-treponemal slide agglutination test for the qualitative and semiquantitative detection of plasma reagins in human serum. Carbon particles coated with a lipid complex are agglutinated when mixed with samples containing reagins.

CLINICAL SIGNIFICANCE

Reagins are a group of antibodies against some components of the damage tissue from patients infected by treponema pallidum, the agent which causes the syphilis. This microorganism produces some damage to the liver and heart, releasing some tissue fragments. Immunological patient system reacts producing reagins, antibodies against these fragments.

REAGENTS

R1 Carbon particles coated with a lipid complex, cardiolipin, lecithin and cholesterol, in phosphate buffer 20 mmol/l, pH 7.0

R2 Positive control red cap: human artificial serum with regain titer $\geq 1/4$

R3 Negative control blue cap: animal serum. Preservative.
- Slide with 6 divisions
- Stirrer

SAMPLES

Fresh serum or plasma. Stable 7 days at 2-8°C or 3 months at -20° C. Samples with presence of fibrin should be centrifuged. Do not use highly hemolyzed or lipemic samples.

PRECAUTION IN USE

 Components from human origin have been tested and found to be negative for the presence of HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

CALIBRATION

The reagent sensitivity is calibrated against the "Human Reactive Serum" from CDC (Centre for Disease Control).

PREPARATION

RPR carbon: swirl the reagent gently to disperse the carbon particles before use. Open the RPR carbon vial, place the micropipette to the dispensing vial and draw by suction the required volume of RPR carbon. Once the test is completed, return the reagent to the original vial and rinse the micropipette and vial with distilled water.

STORAGE AND STABILITY

All the kit components will remain stable until the expiration date printed on the label, when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not freeze: frozen reagents could change the functionality of the test

Always keep vials in vertical position. If position is changed, gently mix to dissolve aggregates that may be present.

Reagents deterioration: presence of particles and turbidity.

ADDITIONAL EQUIPMENT

Mechanical rotator with adjustable speed at 80/100 r.p.m.
Vortex mixer and micropipette from 20/50 µl.

QUALITATIVE PROCEDURE

1. Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50µl of the sample and one drop of each positive and negative controls into separate circles on the slide test.
3. Mix the RPR carbon reagent vigorously or on a vortex mixer before using. Invert the dropper assembly and press gently to remove air bubbles from the micropipette.
4. Place the micropipette in a vertical position and perpendicular to the slide, and add one drop (20 µl) of this reagent next to the sample to be tested.
5. Mix the drop with a stirrer, spreading them over the entire surface of the circle. Use different stirrer for each sample.
6. Place the slide on a mechanical rotator at 80-100 r.p.m. for 8 minutes. False positive results could appear if the test is read later than two minutes.

SEMIQUANTITATIVE PROCEDURE

1. Make serial two fold dilutions of the sample in 9 g/l saline solution.
2. Proceed for each dilution as in the quantitative method.

READING AND INTERPRETATION

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator. Rotate the slide twice by hand before reading.

Interpretation

| Agglutination | Reading | Report |
|--|---------|-----------------|
| Medium or large clumps | R | Reactive |
| Small clumps | W | Weakly reactive |
| No clumping or very slight "roughness" | N | Non reactive |

The titer, in the semiquantitative method, is defined as the highest dilution showing a positive result.

QUALITY CONTROL

Positive and negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation.

All results different from the negative control will be considered as a positive.

PERFORMANCE CHARACTERISTICS

1. Analytical sensitivity: Accurate titer determination of the reference material, under the described assay conditions (see "Calibration").
2. Prozone effect: no prozone effect was detected up to titer $\geq 1/128$.
3. Diagnostic sensitivity: 100%
4. Diagnostic specificity: 98%

INTERFERENCES

Bilirubin does not interfere up to concentration of 20 mg/dl.

Hemoglobin does not interfere up to concentration of 10 g/l.

Lipemia does not interfere up to concentration of 10 g/l.

Rheumatoid factors interfere up to concentration of 300 IU/ml.

LIMITATIONS OF THE PROCEDURE

- High temperature may cause test components to dry on the slide giving an agglutination aspect that can be interpreted as false positive results. It is recommended to place the slide under a humidifying cover.

- RPR carbon test is non-specific for syphilis. All reactive samples should be retested with treponemal methods such as TPHA and FTA-Abs to confirm the results.

- A non reactive result by itself does not exclude a diagnosis of syphilis. Clinical diagnosis should not be made by findings of a single test result, but should integrate both clinical and laboratory data.

- False positive results have been reported in diseases such as infectious mononucleosis, viral pneumonia, toxoplasmosis, pregnancy and autoimmune diseases.

BIBLIOGRAPHY

Gorge P.Schimid. Current Opinion in Infectious Diseases 1994;7 34-40.

Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.

Sandra Larsen et al. A manual of test for Syphilis American Public Health Association 190: 1-192.

Joseph Earle Moore et al. Gastrointestinal Haemorrhage 1952 ; 150 (5) : 467-473.

Young DS. Effects of drugs on clinical laboratory test, 4th ed AACCC Press. 1995.

SYMBOLS

 Red instruction for use  Biological risk

 CE mark (requirement of 98/78 CE regulation)

 Storing temperature limits

 In vitro medical device  Producer