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Fabbricante

Rapid test on card for HCG detection in serum, plasma and urine

PRINCIPLE

Human Chorionic gonadotropin (HCG) is a glycoprotein hormone secreted by the developing placenta shortly after fertilization. In normal pregnancy HCG may become detectable in serum since 7 days after conception and double every 1, 2-3 days. At the first missed menses the HCG concentration is about 100 mIU/ml, but the maximum value of 100-200.000 mIU/ml is detectable at the end of third month. After seven days until 100-200.000 mIU/ml at the end of the third month. HCG is a rapid test for qualitative detection in serum and urine. The method, using a combination of monoclonal antibodies conjugated with cromogen and polyclonal antibodies on membrane, allows to detect HCG levels greater than 25 mIU/ml within 5 minutes.

REAGENTS

Cards
Disposable plastic dropper

STORAGE AND STABILITY

- Cards are stable up to expiry date indicated on the label, stored at room temperature (15-25°C) or refrigerated (4-8°C).
- Do not freeze the test kit.

PREPARATION OF REAGENTS

Reagents are ready to use.

SAMPLE

Serum, plasma and urine.

Serum:

- Collect whole blood into collection tube containing no anticoagulants under the standard laboratory conditions (aseptically in such a way as to avoid haemolysis).
- Allow the blood to clot at room temperature about 30 minutes and separate the serum by centrifugation for 5-10 minutes (at 3000rpm).
- Carefully withdraw the serum for testing, or label and store at 4-8°C for up to seventy two hours.

Urine:

- The urine specimens collected at any time of day may be used; however, the first morning urine generally contains the highest concentration of HCG.
- Collect the urine specimen in a clean container without preservatives.
- If specimens are not to be tested they should be refrigerated immediately at 4 to 8°C. If storage periods greater than 3 days the specimen should be frozen. If samples are refrigerated, they must be equilibrated to room temperature before testing. Avoid repeated freezing and thawing.
- Urine samples exhibiting visible precipitates should be filtered, centrifuged or allowed to settle and clear aliquots are obtained for testing.

Plasma:

- Collect whole blood into a purple, blue or green top collection tube (containing EDTA, citrate or heparin, respectively).
- Separate the plasma by centrifugations.
- Carefully withdraw the plasma for testing, or label and store at 4-8°C for up to seventy two hours.


PROCEDURE

- Bring cards and specimens at room temperature before testing.
- Do not open pouches until ready to perform the assay.
- Using the dropper provided, dispense 2-3 drops (~90-135 µl) of serum, plasma or urine in the sample well (S).
- Wait for red-colored bands to appear. Depending on the concentration of HCG, positive results may be observed as soon as 60 seconds. However, low values of HCG may require more time.


IMPORTANT To confirm negative results with patient's samples the complete reaction line of 5 minutes is required. Do not interpret results after 10 minutes.

READING AND INTERPRETRATION

POSITIVE: Two red lines (bands) are visible in the control ("C") and test ("T") areas of the test window. The intensity of the test line may be less than that of the control line; this still means positive result.


NEGATIVE: The control line appears in the test window, but the test line is not visible.

ERROR:

The test is invalid if the control line is not visible at five minutes. The test failed, or the test procedure was not followed properly. Verify the test procedure and repeat the test with a new testing device.

NOTE

HCG level of 25 mIU/ml can be reached shortly after 7-10 days since conception. The maximum HCG level is about after 8 weeks since last menses and decrease during pregnancy. An high HCG level (> 200.000 mIU/ml) can give a wrong negative result. In case of suspicious pregnancy it's advisable repeat the test diluting the sample 1:1. This is a qualitative test. The shade of red in the test band region (T) will vary depending on the concentration of HCG present. However, neither the quantitative value nor the rate of in HCG can be determined by a qualitative test.

PERFORMANCES

Accuracy Correlation

Correlation with Qualitative Visual Test: 102 randomly selected urine, serum or plasma samples were analyzed in parallel with a commercially available qualitative visual HCG test. The results indicated complete agreement (57 positive samples and 45 negative samples).

102 samples (urine, plasma, serum)	Tests	
	SGM	Qualitative Visual Test
+	57	57
-	45	45

Additionally, 20 urine and 20 serum samples at the HCG concentrations of 0,25,40,80,100mIU/ml were blind labeled and tested in other medical laboratory. The results from this study gave 100% agreement with the expected results.

Sensitivity

- Test sensitivity is 25 mIU/ml according to WHO 3rd International Standard. Serum, plasma or urine specimens containing less than 25mIU/ml HCG may also produce a faint positive result.

To evaluate the sensitivity of the Test at low levels of HCG, the following experiments were carried out: Urine, serum or plasma samples from 100 known non-pregnant subjects were spiked with HCG to the concentrations of 0, 10, 25,40, 100mIU/ml. A total of 100

	HCG concentration (mIU/ml)				
	0	10	25	40	100
Samples	20	20	20	20	20
Positive	0	10	20	20	20
Negative	20	10	0	0	0

samples were blind labeled and tested with the HCG test. Results are summarized in the table below.

- Urine from healthy and non-pregnant women contain HCG levels not detectable by this test.

- This test is capable to detect pregnancy by the first day of the missed menstrual period, although the HCG levels are variable in the early pregnant women.

Specificity

Cross reactivity

The specificity of the HCG Test was determined from cross-reactivity studies with known amounts of homologous hormones.

The following concentration of homologous hormones are found to have no interferences with HCG test:

hTSH	1000 uIU/ml
hLH	300 mIU/ml
hFSH	1000 mIU/ml

Interferences

Potentially interfering substances were added to HCG-free and 25mIU/ml HCG-spiked urine or serum samples. At the concentrations tested, none of the substances interfered in the assay.

Acetaminophen	20 mg/dl
Acetylsalicylic acid	20 mg/dl
Ascorbic acid	20 mg/dl
Atropine	20 mg/dl
Caffeine	20 mg/dl
Centesic acid	20 mg/dl
Glucose	2 g/dl
Haemoglobin	20 mg/dl

LIMITATIONS OF THE TEST

- In addition to pregnancy, HCG has been found in patients with trophoblastic disease. Since the HCG of trophoblastic neoplasms is similar to that found in pregnancy, these conditions, which include

choriocarcinoma and hydatichiform mole, should be ruled out before a diagnosis of pregnancy is reached.

- A normal pregnancy can be distinguished from an ectopic pregnancy based on HCG levels alone. Also, spontaneous miscarriage may cause confusion in interpreting test results.

- A very early pregnancy containing an extremely low concentration of HCG can give a negative result. In this case, repeat the test after at least 48 hours.

- HCG levels remain detectable for several weeks after normal delivery, delivery by caesarean section, spontaneous or therapeutic abortion.

- Some serum specimens with high rheumatoid factor concentration (RF) may yield non specific positive results during testing. Such cases should be discriminated before testing.

- Immunologically interfering substances such as those used in antibody therapy treatments may invalidate the test result.

- The physician should evaluate data obtained by the use of this kit in light of other clinical information.

WASTE MANAGEMENT

Please refer to the local legal requirements.

Note

The number indicated on the kit components identify the production lot. The number indicated on the external label identify the kit package. Any complaint must be furnished with one of these numbers.

BIBLIOGRAPHY

- Batzer, F.R. *Fertility and Sterility*, Vol.34,1,1980.
 Catt, K.J., Dufan, M.L. and Vaitukaitis, J.L.J. *Clin. Endocrinol. Metab.*, Vol.40,537,1975.
 Baunstein, G.D., Rasor, J., Adler D., Danzer, H., Wade, M.E. *Am.J. Obstet. Gynecol.*, Vol.126,678,1976.
 4. Lenton, E.A., Neal, L.M., Sulaiman, R. *Fertility and Sterility*, Vol.37,773,1982.
 Willcox, A.J. Weinberg, C.R., O'Connor J.F., Baird D., D Schlatterer, J.P. Canfield R.E., Armstrong E.G. Nisula, B.C. *Incident of Early Loss of Pregnancy. N.Eng.J. Med.*, Vol.319,189-194,1988.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer