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BIBLIOGRAFIA

Baselt, R.C. Disposition of toxic Drugs and Chemical in Man, Davis, Biomedical Communication, 1982.
 Urine teting for Drugs of Abuse, National Institute on Drugs Abuse (NIDA), Research Monograph 73, 1986.
 Fed. Register, Department of Health and Human Services, mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988.
 McBay, A.J. Clin. Chem. 33, 33B-40B, 1987.
 The Pharmacological Basis of Therapeutics, Gilman, A.G. Goodman, L.S. and Gilman A. eds. MacMillan Publishing, New York, NY, 1980.

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Fabbricante

Rapid test on card for amphetamine detection in urine

PRINCIPLE

Amphetamine is a potent substance also used for therapeutic purpose, with sympathetic-mimetic characteristics. Heavy abuse cause stimulation of central nervous system giving euphoria, alert and increased energy. Amphetamine affect cardiovascular system giving increase blood pressure, cardiac arrhythmias and also anxiety, paranoia and psychotic behaviour. Amphetamine is eliminated in urine in hydroxilated form as deaminated derivatives or integral molecule. The test use an immunometric "competitive" method for amphetamine detection by competition between amphetamine eventually present in the tested sample and amphetamine conjugate fixed to the membrane for a limited number of antibodies anti-amphetamine present in the conjugate anti-amphetamine-colouring. A sufficient quantity of amphetamine, above cut-off levels, will occupy all binding sites of antibodies and will prevent the bond between the conjugate anti-amphetamine-colouring and the coloured band formation on the membrane. Therefore, the formation of one coloured band shows a positive assay. Instead, the contemporaneous appearance of two coloured bands indicates a negative result. This method represents a valid, simple and fast alternative to classic methods as gas chromatography (GC) or mass spectrometry (MS) which require complicated equipment and longer execution times. However, this assay provides only a preliminary analytical test result and, in case of positivity, it must be confirmed by other chemical researches.

REAGENTS

Cards (AMP)	10
Specimens pipettes	10

STORAGE AND STABILITY

Cards are stable up to expiry date showed on the package, stored at room temperature (15-25°C) or refrigerated (4-8°C) .

PREPARATION OF REAGENTS

Reagents are ready to use.

SAMPLE

Fresh urine.

QUALITY CONTROLS

Good laboratory practice recommends the use of control materials to ensure proper kit performances. When testing the positive and negative controls, use the same assay procedure as with a urine specimen.

PROCEDURE


- Allow specimens and cards to come to room temperature prior to testing.
- Fill the urine dropper with specimens and, by holding it vertically, dispense 6 drops (200 µl) of urine in the sample well A of the card.
- Wait 3-5 minutes to read the results. Do not read after 8 minutes because the intensity of the coloured bands may change, or a new line may appear which was no present at the end of the incubation time (5 minutes). Therefore, it's most important to record the results of assay immediately after the five minutes incubation period.

Prudences


- Urine specimens and all materials coming in contact with them should be handled and disposed as if it is infectious and capable of transmitting infection. Avoid contact with skin.
- Avoid cross-contamination of urine samples by using a new specimen cup and specimen pipette for each urine sample.
- Test the sample as soon as possible, preferably the same day of specimen collection.
- The urines may be refrigerated at 2-8°C for 2 days or frozen at -20°C for longer periods.
- Specimen that were refrigerated must reach room temperature (15-25°C) prior to testing.
- Specimen previously frozen must be thawed and mixed thoroughly prior to testing.

INTERPRETATION OF RESULTS

POSITIVE: One colour band appears on the control window (C). The amphetamine concentration is over cut-off levels.

NEGATIVE: Two colour bands appear, even if with weak intensity. The amphetamine concentration is under cut-off levels.

ERROR: If there is no distinct colour band visible both in the test window and test control window the test is inconclusive. It is recommended to repeat the test.

Note

- A very faint line on the test region indicates that the amphetamine in the sample is near or below the cut-off level for the test. However, if any line is observed at the test region, the sample is negative.
- Confirm eventually positive results by opportune confirmation assay.

ANALYTICAL PERFORMANCES

Specificity

- The test is specific for amphetamine and metylendioxyamphetamine.

Interferences

- The following compound was found to cross-react with Amphetamine:

Compound	Concentration (ng/ml)
D-Amphetamine	1,000
1-Amphetamine	10,000
(+/-)3,4-methylenedioxyamphetamine (MDA)	5,000

- The following compounds were found not to cross-react when tested at concentrations up to 100 µg/ml:

Acetone	Acetaminophen
Albumin	Amitriptyline
Ampicillin	Aspartame
Aspirin	Atropine
Benzocaine	Bilirubin
Caffeine	Chloroquine
(+)-Chlorpheniramine	(+ -)Chlorpheniramine
Creatinine pheniramine	Dextromethorphan
Dexbrom	4-Dimethylaminoantipyrine
Doamine	Tyramine
(+/-)-Ephedrine	(-)-Epinephrine
Erythromycin	Ethanol
Furosemide	Guaiacol Glyceril Ether
Glucose	Hemoglobin
Imipramine	(+ -)-Isoproterenol
(1R,2S)-(-)-N-Methyl-Ephedrine	Lidocaine
d-Methamphetamine	1-Methamphetamine
(+/-)3,4-methylenedioxymethamphetamine	
(+)-Naproxen	(+/-)-Norephedrine
Oxaloc Acid	Penicillin G
Pheniramine	Phenothiazine
1-Phenylephrine	B-Phenylethylamine
Procaine	Quinidine
Ranitidine	Riboflavin
Sodium Chloride	Sulindac
Thioridazine	Trifluoperazine
Trimethobenzamide	Tyramine
Vitamin C	

Correlation-accuracy

The comparison between the results obtained by the kit amphetamine easy-card and a commercial immunoenzymatic kit on 55 positive samples and 92 negative samples, gives a superposition of results at 99%. The correlation on 55 samples, resulted positives by GC/MS assays, gives a superposition of results at > 99%.

Sensitivity

The test sensitivity is ≥ 1000 ng/ml.

LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates only the presence of a drug/metabolite and does not indicate or measure intoxication.
- There is a possibility that other substances and/or factors not listed may interfere with the test and cause false results, such as technical or procedural errors. See "Specificity" for lists of substances that will produce positive results or that do not interfere with tests performance.

WASTE MANAGEMENT

Please refer to the local legal requirements.

Note

The number showed on single product refer to the production lot. The number showed on the exterior label refer to the packaging. In case of claims, please refer to one of these two number.

BIBLIOGRAPHY

Baselt, R.C. Disposition of toxic Drugs and Chemical in Man, Davis, Biomedical Communication, 1982.
Urine teting for Drugs of Abuse, National Institute on Drugs Abuse (NIDA), Research Monograph 73, 1986.
Fed. Register, Department of Health and Human Services, mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988.
McBay, A.J. Clin. Chem. 33, 33B-40B, 1987.
The Pharmacological Basis of Therapeutics, Gilman, A.G. Goodman, L.S. and Gilman A. eds. MacMillan Publishing, New York, NY, 1980.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer