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### BIBLIOGRAFIA

Baselt, R.C. Disposition of toxic Drugs and Chemical in Man, Davis, Biomedical Communication, 1982.  
 Urine teting for Drugs of Abuse, National Institute on Drugs Abuse (NIDA), Research Monograph 73, 1986.  
 Fed. Registerp, Department of Health and Human Services, mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988.  
 McBay, A.J. Clin. Chem. 33, 33B-40B, 1987.  
 The Pharmacological Basis of Therapeutics, Gilman, A.G. Goodman, L.S. and Gilman A. eds. MacMillan Publishing, New York, NY, 1980.

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Fabbricante

## Rapid test on card for opiates detection in urine

### PRINCIPLE

Opiates analgesics include a big group of compounds which act as depressives on central nervous system. Morphine characterizes these compounds and it is eliminated not metabolised after the assumption of codeine and heroin. Morphine is present in urine also many days after the assumption of an opiates dose.

The test uses a "competitive" immunometric method for morphine detection by the competition between morphine fixed on membrane and morphine free, eventually present in the tested sample, to the conjugated anti-morphine-colouring. A sufficient quantity of morphine, above cut-off level, will occupy all binding sites of antibodies and will prevent the bond between conjugated anti-morphine-colouring and the coloured band formation on the membrane. Therefore, the formation of one coloured band shows a positive assay. Instead, the contemporaneous appearance of two coloured bands indicates a negative result. This method represents a valid, simple and fast alternative to classic methods as gas chromatography (GC) or mass spectrometry (MS) which require complicated equipment and longer execution times. However, this assay provides only a preliminary analytical test result and, in case of positivity, it must be confirmed by other chemical researches.

### REAGENTS

Cards (MOP)	10
Specimens pipettes	10

### STORAGE AND STABILITY

Cards are stable up to expiry date showed on the package, stored at room temperature (15-25°C) or refrigerated (4-8°C).

### PREPARATION OF REAGENTS

Reagents are ready to use.

### SAMPLE

Fresh urine.

### QUALITY CONTROLS

Good laboratory practice recommends the use of control materials to ensure proper kit performances. When testing the positive and negative controls, use the same assay procedure as with a urine specimen.

### PROCEDURE


- Allow specimens and cards to come to room temperature prior to testing.
- Fill the urine dropper with specimens and, by holding it vertically, dispense 6 drops (200 µl) of urine in the sample well A of the card.
- Wait 3-5 minutes to read the results. Do not read after 8 minutes because the intensity of the coloured bands may change, or a new line may appear which was not present at the end of the incubation time (5 minutes). Therefore, it's most important to record the results of assay immediately after the five minutes incubation period.

### Prudences


- Urine specimens and all materials coming in contact with them should be handled and disposed as if it is infectious and capable of transmitting infection. Avoid contact with skin.
- Avoid cross-contamination of urine samples by using a new specimen cup and specimen pipette for each urine sample.
- Test the sample as soon as possible, preferably the same day of specimen collection.
- The urines may be refrigerated at 2-8°C for 2 days or frozen at -20°C for longer periods.
- Specimens that were refrigerated must reach room temperature (15-25°C) prior to testing.
- Specimens previously frozen must be thawed and mixed thoroughly prior to testing.

### INTERPRETATION OF RESULTS

**POSITIVE:** One colour band appears on the control window (C). The morphine concentration is over cut-off levels.

**NEGATIVE:** Two colour bands appear, even if with weak intensity. The morphine concentration is under cut-off levels.

**ERROR:** If there is no distinct colour band visible both in the test window and test control window the test is inconclusive. It is recommended to repeat the test.

### Note

- A very faint line on the test region indicates that the morphine in the sample is near or below the cut-off level for the test. However, if any line is observed at the test region, the sample is negative.
- Confirm eventually positive results by opportune confirmation assay.

### ANALYTICAL PERFORMANCES

#### Specificity

- The test is specific for morphine and opiates in urine.
- The test specificity was tested by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

### Interferences

- The following compounds produce positive results when tested at levels equal to or greater than the concentrations listed below.

Compound	Concentration (ng/ml)
Codeine	300
Ethyl morphine	300
Hydromorphone	400
Morphine-3-β-d-glucuronide	490
Morphine	300
Hydrocodone	375
Thebaine	30,000

- The following compounds were found not to cross-react when tested at concentrations up to 100 µg/ml:

Acetaminophen	Imipramine
Cyclobenzaprine	Phencyclidine
Acetone	(+ -)-Isoproterenol
r-Cyclodextrin	Phenelzine
Albumin	Lidocaine
C-yproheptadine	Pheniramine
Amitriptyline	Maprotiline
(-)-Desoxyephedrine	Phenobarbital
D-Amphetamine	Mephentermine
Dextromethorphan	Penothiazine
L-Amphetamine	MDA
Diazepam	Meperidine
Amaxaoine	Mephentermine
4-Dimethylaminoantipyrine	Ampicillin
β-Phenylethylamine	Diphenhydramine
Mes	Aspartame
Phentermine	(Asp-Pheethyl-Ester)
Methadol	5,5-Diphenylhydantoin.
Primidone	Aspirin
Methadone	(Acetylsalicylic Acid)
Procaine	Dopamine
Dopamine	Methamphetamine
(3-Hydroxytyramine)	Promazine
Atropine	Methapyrilene
Doxylamine	Promethazine
Baclofen	Methaqualone
Ecgonine	d-Propoxyphene
(1R,2S)-(-)-N-Methyl-ephedrine	Ecgonine Methy Ester
Benzocaine	Protriptyline
(Ethyl-Aminobenzoate)	Methylphenidate
(-) Ephedrine	Pseudoephedrine HCL
Benzoyllecgonine	(+/-)3,4-Methylenedioxy-
(+/-)-Ephedrine	2-propylantanoic acid
Bilirubin	Methamphetamine(+/-)3,4-MDMA
(-)- Ephedrine	Quinidine
(+)-Brompheniramine	Naloxone
(Dexbrompheniramine)	Quinine
(+)- Ephedrine	Naltrexone
Caffeine	Ranitidine
β-Naphthaleneacetic acid	Carbamate
Riboflavin	Erythromycin
(+)-Naproxen	Carbamazepine
Salicylic acid	Ethanol
(-)-Nicotine	Carisoprodol
Secobarbital	Chloroquine

Nicotinic acid	(+ -)Chlorpheniramine
Sodium Chloride	Fentanyl
(+/-)-Norephedrine	(+/-)-Chlorpheniramine
((+/-)-Phenylpropanol)	Fenofibrate
Chlorprothixene	Nortriptyline
Gemfibrozil	Theophylline
Chlorthalidone	Noscapine
Glucose	cis-Thiothixene
Cocaine	Oxalic acid
Guaiacol glyceryl ether	Trifluoperazine
Clofibrate	Oxycodone
DL-Homatropine	Trimethobenzamide
Creatine	Penicillin G
Hemoglobin	Tyramine
Creatinine	Pentobarbital
Ibuprofen	Vitamin C
Perphenazine	

### Correlation-accuracy

The comparison between the results obtained by the kit opiates easy-card and a commercial immunoenzymatic kit on 62 positive samples and 178 negative samples, gives a superposition of results at 100%. The correlation on 73 samples, resulted positives by GC/MS assays, gives a superposition of results at > 99%.

### Sensitivity

The test sensitivity is  $\geq 300$  ng/ml.

### LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates only the presence of a drug/metabolite and does not indicate or measure intoxication.
- There is a possibility that other substances and/or factors not listed may interfere with the test and cause false results, such as technical or procedural errors. See "Specificity" for lists of substances that will produce positive results or that do not interfere with tests performance.

### WASTE MANAGEMENT

Please refer to the local legal requirements.

### Note

The number showed on single product refer to the production lot. The number showed on the exterior label refer to the packaging. In case of claims, please refer to one of these two number.

### BIBLIOGRAPHY

Baselt, R.C. Disposition of toxic Drugs and Chemical in Man, Davis, Biomedical Communication, 1982.  
 Urine testing for Drugs of Abuse, National Institute on Drug Abuse (NIDA), Research Monograph 73, 1986.  
 Fed. Register, Department of Health and Human Services, mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988.  
 McBay, A.J. Clin. Chem. 33, 33B-40B, 1987.  
 The Pharmacological Basis of Therapeutics, Gilman, A.G. Goodman, L.S. and Gilman A. eds. MacMillan Publishing, New York, NY, 1980.

### SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer