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Fabbricante

Rapid test on card for antibodies anti-Helicobacter Pylori detection in serum, plasma and blood

PRINCIPLE

Helicobacter Pylori test allows to detect qualitative IgG anti-H. Pylori in serum, plasma or whole blood. H. Pylori is a gram negative bacteria, characterized completely only in the last ten years by cellular cultures. This bacteria can colonize gastric mucous wall for long periods, even several years, and their presence is linked with chronic gastritis, superficial or diffused. The epidemiologic characteristics of gastric carcinoma are similar to those coming from H. Pylori infections and moreover recent researches have highlighted that gastritis from H. Pylori are a risk factor for gastric carcinoma. The classical analysis methods for direct research of H.Pylori are based on samples collection by gastric endoscopy and on microorganisms research by technique of cellular cultures or histological. At the moment, the method for the research of H. Pylori infection is based on serologic research of anti- H. Pylori (these methods are easy to use and allow to follow the infection in the different stages in patients under antibacterial therapy) because there is an high correlation, in infectious patients, between the presence of these antibodies in serum and in the stomach. The test is a chromatographic immunoassay for the qualitative determination of anti-H.Pylori IgG antibodies in human serum, plasma or whole blood. The test device contains a membrane strip which is pre-coated with H. Pylori antigens on the test band region and specific monoclonal antibody on the control band region. The H.P. antigens-colloid gold conjugate pad is placed at the end of the membrane. When IgG anti-H.P. are present in patient sample, the mixture of colloid gold conjugate, patient sample and buffer moves along the membrane to the test region T and form a visible line. This indicates a positive result for the detection of specific IgG antibodies. When IgG anti-H.P. are absent in the sample, no visible band will appear in the T region. This indicates a negative result for the detection of H.P. specific IgG antibodies. A line will always appear in the C region showing that the device is correctly working.

REAGENTS

R1	Cards	10
R2	Development solution	2 ml

STORAGE AND STABILITY

- Cards are stable up to expiry date showed on the package, stored at room temperature (15-25°C) and in the sealed pouch.
- Do not freeze the kit.

PREPARATION OF REAGENTS

Reagents are ready to use.

SAMPLE

Serum, plasma and whole blood (EDTA, citrate o heparin)
Do not use whole blood which have been stored for more than three days. Serum can be stored at 4-8°C for up to two weeks and frozen at -20°C up to one year. Plasma can be stored at 4-8°C for up to one year. Heated sample may cause erroneous results.

WASTE MANAGEMENT


Please refer to the local legal requirements.

PROCEDURE


- Allow specimens and cards to come to room temperature prior to testing.
- Fill the dropper with specimens and, by holding it vertically, dispense 1 drop (25 µl) of serum and plasma or 2 drops (40µl) of whole blood in well S of the card.
- Add 2-3 drops of development solution in the sample well S.
- Read the results at 5 minutes Do not interpret the results after 8 minutes. Be sure that the background of the card is clear.

INTERPRETATION OF RESULTS

POSITIVE: in presence of a colored band in the well **T** confirmed by a colored band in the control well **C**.

NEGATIVE: in absence of a colored band in well **T** and the presence of a colored band in the well control **C**.

INVALID: If there is no distinct colour band visible both in the test window and test control window the test is inconclusive. It is recommended to repeat the test.

QUALITY CONTROL

A procedural control is included in the test. A colored band will appear on the control region C showing proper performances of the device.

ANALYTICAL PERFORMANCES

Accuracy – comparison with biopsy results

The accuracy of the kit was evaluated in comparison to biopsy results of human specimens over 317 samples. This comparison study results gave a sensitivity of 92.6% (100/215), a specificity of 82.4% (84/102), and a total agreement of 89.3% (283/317). The relatively low specificity of the serological test results in comparison to the biopsy results may be partially attributed to a sampling error of the biopsy test.

Comparison study with a commercially available H.P. test

The accuracy of this kit was evaluated against a commercially available rapid test using 170 serum/plasma specimens. These results gave a sensitivity of 84.0% (63/75) a relative specificity of 100% (95/95) and a total agreement of 92.9 % (158/170).

Specificity

The cross-reactivity of the kit was evaluated by an inhibitory assay. Testing was performed with Campylobacter coli, fetus, jejuni and E. Coli. Different antigens were added to patient samples and tested. All species tested showed no cross reactivity. Based on these results, it can be concluded that this device is specific for H.Pylori antibodies.

LIMITATIONS OF PROCEDURE

This kit is to be used for the qualitative detection of IgG antibodies to H.Pylori. The kit should be used for symptomatic individuals suspected of having gastrointestinal disorder. Diagnosis of gastritis and/or peptic ulcer should be made by confirmation with other clinical findings. A positive results suggest the presence of IgG antibodies to H.P.; it does not distinguish between active infection and past exposure to H.Pylori and does not necessarily indicate the presence of gastrointestinal disease. A negative result does not rule out H. Pylori infection because the IgG antibodies may not be present in sufficient quantity to be detected.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer