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BIBLIOGRAFIA

Barbee R. A., mM. Halonen, M. Lebowitz, and B. Burrows. 1981. Distribution of IgE in a community population sample; correlations with age, sex, and allergen skin reactivity. *J. Allergen Clin. Immunol.* 68: 106-111.

Kjellman, N-I.M, S.G.O. Johansson, and A. Roth. 1976. Serum IgE levels in healthy children quantified by a sandwich technique (PRIST). *Clin. Allergy* 6: 51-59.

Ishika, T. (1982), IgE and Mechanisms of Ige-mediated hypersensitivity, *Ann. Allergy* 46, 313-319.

Geha, RS (1984), Human IgE, *J Allergy & Clin Immunol*, 74, 109-120.

Witting, HJ, Belliot, J, Dfe Fillipi; I and Royal, G (1980), Age-related serum Immunoglobulin E levels in healthy subjects and in patients with allergic disease, *J Allergy & Clin Immunol*, 66, 305-313.

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visual card

For in vitro medical device
31039 - 10 test

Rapid test on card for IgE detection in serum

PRINCIPLE

Immunoglobulin E (IgE) is an immunoglobulin normally found in nanogram level in healthy individual. In allergic patients the IgE level is increased.

Screening for the total IgE is of diagnostic use in case of asthma (to differentiate allergic from non-allergic asthma), rhinitis and eczema. Other disease in which elevated serum IgE levels are found include parasitic infections, broncho-pulmonary aspergillosis and some of dermatitis. The IgE card is a rapid, qualitative assay for the detection of IgE in serum. The method employ a "capture" immunoassay on membrane for IgE detection. In the membrane there is an antigen anti-IgE conjugated with a colouring agent, this one binds to the IgE forming an antibody-antigen complex. This complex migrates on the membrane, and the test positivity is showed by a rose pink band appearance. The reaction mixture continues through the control zone where should be appeared a rose pink band demonstrating that the reagents are functioning correctly. The test is sensitive to an IgE concentration of 150 IU/ml.

REAGENTS

Cards (IgE) 10
Specimens pipettes 10

STORAGE AND STABILITY

- Reagents, if stored at room temperature (15-25°C) or refrigerated (4-8°C) in sealed package, are stable up to expiry date indicated.
- Do not use reagents over the expiry date indicated on the package.
- Do not freeze the kit.

PREPARATION OF REAGENTS

Reagents are ready to use.

SAMPLE

Serum.

LIMITI DEL METODO

- The test is specifically designed to detect IgE in serum sample.
- Data obtained should be evaluated in light of other clinical information.

WASTE MANAGEMENT

Please refer to the local legal requirements.

Note

The number showed on single product refer to the production lot. The number showed on the exterior label refer to the packaging. In case of claims, please refer to one of these two number.

BIBLIOGRAPHY

Barbee R. A., mM. Halonen, M. Lebowitz, and B. Burrows. 1981. Distribution of IgE in a community population sample; correlations with age, sex, and allergen skin reactivity. J. Allergen Clin. Immunol. 68: 106-111.

Kjellman, N-I.M, S.G.O. Johansson, and A. Roth. 1976. Serum IgE levels in healthy children quantified by a sandwich technique (PRIST). Clin. Allergy 6: 51-59.

Ishika, T. (1982), IgE and Mechanisms of Ige-mediated hypersensitivity, Ann. Allergy 46, 313-319.

Geha, RS (1984), Human IgE, J Allergy & Clin Immunol, 74, 109-120.

Witting, HJ, Belliot, J, Dfe Fillipi, I and Royal, G (1980), Age-related serum Immuyoglobulin E levels in healthy subjects and in patients with allergic disease, J Allergy & Clin Immunol, 66, 305-313.

PROCEDURE


- Allow specimens and cards to come to room temperature prior to testing.
- Fill the dropper with sample and, by holding it vertically, dispense 5 drops (200 µl) in the sample well A of the card.
- Read the result after 5-10 minutes.

Prudences


- The serum specimen should be collected under the standard laboratory conditions (aseptically in such a way a sto avoid haemolysis)
- No centrifugation or filtration is required
- The specimen can be stored at 4-8°C within 48 hours.
- For longer times, the specimen should be frozen. In this case, the specimen must be thawed, thoroughly mixed and brought to room temperature prior to testing. Avoid repeated freezing and thawing.
- In the case of cloudiness, high viscosity or presence of particulate matter into the serum specimen, it should be diluted 1:1 with diluting buffer before testing.

INTERPRETATION OF RESULTS

POSITIVE: Two coloured bands (B+C). The IgE concentration is over 150 IU/ml.

NEGATIVE: One coloured band appears (C).

ERROR: If there is no distinct colour band visible both in the test window and test control window the test is inconclusive. It is recommended to repeat the test.

ANALYTICAL PERFORMANCES

Sensitivity

The test sensitivity is ≥ 150 IU/ml. Concentration below this value may also be detected occasionally.

Specificity

A study was performed on 55 sera pretested using UNICAP TOTAL IgE (Pharmacia & John). The results are summarized in the following table:

		IgE - CARD	
		+	-
UNICAP	+	31	2
	-	2	20

The 4 discordant results observed between the two methods are due to serum sample with IgE concentration between 100-150 IU/ml, closed to the card cut-off.

From the above table, the test specificity is 90,9% and the correlation is 92,7%.

Hook effect

There was no observed hook effect up to an IgE concentration of 32.000 IU/ml.

Interferences

No reaction was observed with critical RF and complement interfering sera.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer