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Fabbricante

## Rapid Test on card for the detection of albumin in human urine

### PRINCIPLE

The Test is based on the principle of competitive immunochemical reaction between the albumin in the urine sample and albumin immobilized on the membrane for the limited antibody sites. The test strip contains a nitrocellulose membrane strip and a colored anti-human albumin monoclonal antibody-colloidal gold conjugate pad. The nitrocellulose membrane is pre-coated with human albumin on the test band region and a secondary antibody on the control band region. The anti-albumin antibody-colloidal gold conjugate pad is placed at the bottom of the membrane strip. During the test, the urine sample is allowed to migrate upward and dehydrating the anti-albumin antibody-colloidal gold conjugate. The mixture then migrates along the membrane chromatographically by the capillary action to the immobilized albumin zone on the test band region. When albumin is absent in the urine, the colored anti-albumin antibody-colloidal gold conjugate and immobilized albumin bind specifically to form a visible line as the antibody complexes with the albumin. When albumin is present in the urine, the urine albumin competes with albumin immobilized on the test band region for the limited antibody sites on the anti-albumin antibody-colloidal gold conjugate. The test band will become less intense with increasing urine albumin concentration. A control band generated by a different antigen/antibody reaction is also present at control band region of the test strip. The control band should always appear regardless the presence of albumin in the urine sample. Since test band intensity decreases with increasing albumin concentration and the control band intensity is independent of the albumin concentration level in the sample. The presence of control band also serves as a build-in control, which demonstrates that the test is performed properly.

### REAGENTS AND PROVIDED MATERIALS

Reagents are ready to use. Test strip contains human albumin immobilized on the membrane and anti-human albumin monoclonal labeled with colloidal gold.  
Cards (20)

### STORAGE AND STABILITY

The kit should be stored at room temperature or refrigerated (4-30°C) in the sealed pouch.

### PRECAUTIONS

For in vitro diagnostic use.  
Do not use it after expiration date

### SAMPLE

Urine  
First morning urine specimen is preferably to be used for test, and the urine does not require any special handling or pretreatment. Three samples are acceptable: 24-hour urine collection and 12-hour urine collection and overnight urine collection. The specimen may be refrigerated at 2-8°C up to 2 days. Specimens that have been refrigerated must be

equilibrated to room temperature, and mixed thoroughly prior to testing.

### PROCEDURE

Bring the kit and the sample to room temperature (15-25°C) prior to testing.

Use the pipette to withdraw urine from the specimen collection container and dispense 3 drops in a vertical position into sample well.


Read result in 5 minutes. Do not interpret result after 10 minutes.

### INTERPRETATION OF RESULTS

The test results are interpreted by comparing the relative intensity of test band (T) and control band (C). The sample can be classified into one of three categories: Negative (  $\leq 15$ mg/l), Cut-off (15-20mg/l), and positive (>20 mg/l)


#### NEGATIVE ( $\leq 15$ mg/l):

The test band is darker than Control band.


#### Cut-off (15-20mg/l):

The intensities of Test band and Control band are approximately equal.


#### Positive ( $\geq 20$ mg/l):

The test band is lighter than Control band.



#### Dose-hook effect ( $\geq 200$ mg/l)

The test band will disappear.

#### Invalid:

No control band appears in the section of the detection window.

## LIMITATIONS OF PROCEDURE

The assay is designed for use with human urine only.

Dose-hook effect (< 200mg/l):The test band will disappear.

There is a possibility that technical or procedural errors may interfere with the test and cause false results. If questionable results are obtained or abnormal level of microalbuminuria is suspected, the test should be repeated using a fresh urine specimen. Several clinical conditions, such as heavy exercise, urinary tract infection, acute febrile illnesses, and heart failure, may transiently increase urinary albumin excretion (5). Screening of microalbuminuria should be postponed in these situations.

## QUALITY CONTROL

Good laboratory practice should be followed including proper specimen collection and handling, and the use of quality control materials to ensure proper assay performance. Quality control specimens can be purchased from commercial sources. When testing the positive and negative controls, the same assay procedure should be followed.

## PERFORMANCE CHARACTERISTICS

### Accuracy

The accuracy of the One-Step MAU Test was evaluated in comparison to a commercially available immunoassay using 128 urine samples. The result indicated complete agreement (72 positive samples and 56 negative samples).

### Precision

The precision of the One-Step MAU Test was determined by testing two levels of the spiked controls with human serum albumin concentrations of 10 and 50mg/l. The control at 10 mg/l should give a negative result and the control at 50 mg/l should give a positive result. Between-run precision was established by testing two levels of controls 25 times each level over a five-day period. Within-run precision was established by testing two levels of controls 50 times each level.

### Between-run precision

Level(mg/l)	No. Tested	Correct Result	%Correct Result
10	125	125	100
50	125	125	100

### Within-run Precision

Level(mg/l)	No. Tested	Correct Result	%Correct Result
10	50	50	100
50	50	50	100

## Interference Study

Potentially interference of common medications and various human proteins, and other substances that are likely to be present in urine was studied. All substances were spiked into negative and positive controls and tested with One-Step MAU Test. No interference was observed at substance concentrations listed below:

Substances	Concentration (mg/l)
Acetaminophen	20
Acetylsalicylic acid	20
Ampicillin	20
Ascorbic acid	20
Atropine	20
Caffeine	20
Gentisic acid	20
Glucose	2000
Hemoglobin	500
Mestranol	13
Procaine	20
Tetracycline hydrochloride	20

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