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IVD  
REF

For in vitro medical device  
31142 - 40 test

## Rapid test on card for occult blood detection in stool – without diet

### PRINCIPLE

Occult blood test allows the qualitative research of haemoglobin in fecal matter, and can be used as screening test for bleeding of lower intestinal tract which can occur in colon-rectal polyps, carcinoma and Crohn disease.

The advantage of this test is that there is no need for restrictive diets - without meat and food containing peroxidase - to run the test because it is specific for human haemoglobin. The specific antibodies are fixed to the membrane. In presence of haemoglobin forms a coloured band in the test area. A second band confirm the correct test working.

### REAGENTS

- Solid stage: membrane sensitized with specific antibodies for Hb containing a second monoclonal antibody bound with colloidal gold.
- Extraction liquid: highly hypotonic solution.

### STORAGE AND STABILITY

Reagents are stable up to expiry date marked on the box if stored at room temperature.

### PREPARATION FO REAGENTS

Reagents are ready to use.


### SAMPLE

Human fecal matter.

- Haemoglobin in stool is stable up to 7 days at 4°C.

### PROCEDURE

- Take only a little bit amount of stool (20-40 mg) on the collection device by the special stick.
- Screw the extraction vial with the sample collected and shake.
- Break off the tip and dispense 3 drops (150 µl) in the sample hole of the card **S**.
- Wait 5 minutes to read the results. Do not read after 8 minutes.


### PRECAUTION IN USE


- Collect fecal sample in a clean and dry container.
- Open the vial of extraction liquid unscrewing the plug and extracting from it the stick attached.
- Dip the stick into three different points of the same stool sample, only the grained extremity, pull out and insert it in extraction liquid.
- Return the sample loaded into its place in the collection device, firmly screw and shut, then shake.

### INTERPRETATION OF RESULTS

**POSITIVE:** two coloured lines in (T) test hole and (C) control of the card.

**NEGATIVE:** only one coloured line in (C) control hole of the card.

**ERROR:** No coloured line. The test must be repeated.

### Note

Any band appeared after 30 minutes is not important. To preserve the card result longer cover by transparent and adhesive tape the reading window; in this way is avoided the colouring retromigration.

### ANALYTICAL PERFORMANCES

#### Sensitivity

The sensibility limit of the test is 0,5 µg of Hb/g stool corresponding to 0,015 µg of Hb/ml of extraction liquid.

#### Specificity

The test do not measure the haemoglobin from first intestinal tract, because it is modified by acids and digestive enzymes.

### Note

The described test must be considered a rapid screening test for fecal occult blood research that is indicative of gastrointestinal tract diseases associated to bleeding, as colon-rectal carcinoma, colon polyps, Crohn disease and ulcer colitis.

Do not use as the only test for diagnosis of colon-rectal cancer and other specific diseases. To have a complete diagnosis use this test result together to others checks.

The components of this IVD are tested between them without check the compatibility with others producers components. Even if is not excluded the chance that this components are used together the same components of others producer, there is not experimental proof of this compatibility.

### WASTE MANAGEMENT

Please refer to the local requirements.

### Note

The number showed on single product refer to the production lot. The number showed on the exterior label refer to the packaging. In case of claims, please refer to one of these two number.

### BIBLIOGRAPHY

- Parkin, D.M., and al., Int. J. cancer, 1988, 41: 184-197.  
Winaver, S.J., and al., Gastroenterol., 1976, 70: 783-789.  
Winaver, S.J., and al., Gastroenterol., 1977, 72: 1150-1155.  
Frommer, D.J., and al., British medical, 1988, 290: 1092-1094.  
Simon, J.B., Gastroenterol., 1985, 88: 820-837

### SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer